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Expression and purification of GST-OSMR FNIII domains from LB medium

A single colony was picked up from the LB agar plate using a micropipette tip and dipped in 2 ml of LB medium containing ampicillin (100 µg/ml). The culture was incubated at 37 °C for 10–12 h in an orbital shaker at 180 rpm. From this culture, 50 µl of inoculum was transferred into 500ml of LB medium containing ampicillin and grown at 37 °C till the bacterial cell density reached an absorbance of 0.6 at OD₆₀₀. When the required OD was observed, cells were induced with IPTG (1mM) and incubated for another 3 h. Cells were then pelleted down at 4000 rpm, 4 °C, for 15 min. Cell pellets were dissolved in a lysis buffer containing a protease inhibitor cocktail. The cell suspension was sonicated, and the lysate obtained was spun down at 4000 rpm for 15 min at 4 °C. The supernatant obtained was transferred into fresh 15ml centrifuge tubes containing 300 µl of glutathione-agarose beads, rotated at 4 °C for 2 h, and then spun down at 3000 rpm for 5 min at 4 °C. The supernatant was discarded, and the beads were washed thoroughly to get rid of unspecific contaminating proteins. GST-tagged fibronectin domains bound to the beads were eluted by adding 2 ml of elution buffer, mixed thoroughly, and then centrifuged at 3000 rpm, 4 °C, for 5 min (this step was repeated three times to remove the bound protein domains completely). Finally, eluents obtained for each domain were pooled up separately in 15ml centrifuge tubes. Then the eluents were dialyzed (10 KDa cutoff dialysis bag) overnight in a dialysis buffer and concentrated by centrifugation in spin columns at 3000 rpm, 4 °C, for 20 min.

Immunoblotting

Purified fibronectin domains were loaded onto 10% SDS-PAGE and were initially run at 80 V until they crossed the stacking gel, then at a constant voltage of 100 V for approximately one and a half hours. After the completion of SDS-PAGE, western transfer was performed at 100 V for approximately 90 min. After confirming the protein transfer with the ponceau stain, the membrane was blocked for 1 h with 5% skimmed milk prepared in 1X TBST. After blocking, the primary antibody (anti-GST antibody) prepared in 5% skimmed milk (1:1000 dilution) was added to the membrane and incubated for 2 h at room temperature. Later, the membrane was washed thrice with TBST (1X) and incubated with anti-rabbit antibody (secondary antibody) prepared in 5% skimmed milk for 1 h at room temperature. Finally, after the completion of secondary antibody incubation, the membrane was washed thrice with TBST (1X) and processed with enhanced chemiluminescence reagent (Amersham) for observing protein bands.

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
Fig. 4 SDS PAGE of Purified FNIII domains and Immunoblotting: **A** Purified and concentrated GST-tagged FNIII domains were loaded on 10% SDS PAGE. Lane1: Protein Marker, Lane2:GST alone, Lane3:GST-FNIIID1, Lane4:GST-FNIIID2, Lane5: GST- FNIIID3, Lane 6:GST-FNIIID4. **B** The purity of the protein domains was confirmed with Immunoblotting using an anti-GST antibody

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« Sponsored and reviewed by ICCS Quality and Standards Committee^[1-3] Title: Tissue Disaggregation Methods for Flow Cytometric Immunophenotyping »

Written by: Bethany Vallangeon, MD and Ameet R Kini, MD, PhD ; Date: 8/6/2019.

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Sponsored and reviewed by ICCS Quality and Standards Committee

Title: Tissue Disaggregation Methods for Flow Cytometric Immunophenotyping

Written by: Bethany Vallangeon, MD and Ameet R Kini, MD, PhD

Date: 8/6/2019

INTRODUCTION

Successful flow cytometry analysis requires a single-cell suspension; therefore, peripheral blood, bone marrow, and body fluid samples are all very suitable sample types. In contrast, tissue samples such as lymph node and extranodal tissue require processing into single-cell suspension before flow cytometric analysis can be performed. It is advisable that the selected method of tissue processing preserves the cell viability and antigenicity as much as possible. In particular, diagnosis of large cell lymphoma by flow cytometric analysis is often difficult; up to 25% of large cell lymphomas may not be detected.¹

The tissue disaggregation method used can adversely impact immunophenotyping, including possible cell loss, decreased viability, non-specific antibody binding, and variable alterations of antigenicity. It is critical to select a method that is best for your assay and your laboratory and to validate that protocol in house.

Our goal is to compare and contrast different tissue disaggregation methods, discuss the factors that impact the results, advantages and disadvantages, and provide a starting point for choosing and validating the optimal method for your laboratory.

METHODS

The disaggregation of tissue into a single-cell suspension is a critical step in the flow cytometric analysis of hematologic malignancy. Lymphomas dissociate variably based on the nature of involved tissue (extracellular matrix); as an example, lymphomas involving extranodal sites such as skin may present a technical issue. Historically, various methodologies of mechanical dissociation have been utilized including scraping, mincing, fine needle aspiration and automated techniques.

Mechanical Manual Disaggregation Techniques

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In general, manual methods involve using a combination of needle, surgical blade, and/or mesh to disaggregate tissue and washing cells with an appropriate media. The cells can then be collected utilizing a disposable pipette and transferred into a processing tube. The resulting sample can be washed with media or PBS/azide and the cell button resuspended in PBS/azide or media and subsequently counted to adjust the cell concentration.

Figure 1 (below) illustrates the procedure for tissue disaggregation using the needle method while Figure 2 illustrates the mesh technique.


Figure 1 (from Reference 2). Needle Manual Method. A: Puncturing tissue with needle and infusing PBS solution, (B) Aspiration of cell suspension with plastic pipette, (C) Collection tube, and (D) Cell suspension recovery comparison (left tube = manual; right tube = automated).

By George Deeb (Emory University) and Andrea Illingworth (Dahl-Chase Diagnostic Services)
Date: 3/30/2017

Figure 2 (from Reference 3). Mesh Manual Method. In mechanical disaggregation, the tissue sample is disaggregated by simply placing it on a wire mesh over a Petri dish and gently pressing down on the tissue with a glass pestle (or with tissue scissors) and washing cells with appropriate media. The cells can be collected by using a disposable pipette and transferring the cells into a processing tube.

The most prominent advantage of using a manual method is in its notably gentler technique²; though extensive membrane damage due to shear forces have been reported to decrease the number of viable or intact cells.⁶

Vortex disaggregation involves placing a fresh tissue sample in RPMI, then vortexing or vigorously shaking the sample for approximately 10 seconds until the RPMI solution becomes cloudy. Larger specimens may be serially sectioned prior to disaggregation. Small samples, such as GI or skin biopsies may be vortexed in toto. Remaining tissue fragments are subsequently processed for routine histology while the single cell suspension is processed as per any other single-cell suspension for flow cytometric analysis. This method is useful for extremely small samples and “dry-tap” bone marrow biopsies. This is also the only method in which histologic and flow cytometric correlation can be performed on the same tissue sample. The vortex method has been reported to yield lower cell counts from bone marrow core biopsies or extremely small and/or fibrotic samples; however, these samples largely provided diagnostic results and correlated with histology findings.⁶ As with essentially any tissue disaggregation technique, samples with significant amounts of necrosis and/or fibrosis remain a challenge.

Mechanical Automated Disaggregation (Commercial/Medimachine⁵)

Automated methods using commercially available disaggregation machines offer a fast (typically <5 min processing time) and simple technique for creating single cell suspensions from tissue or

bone marrow samples. They require limited tissue handling and provide a closed disposable system with standardized, operator-independent procedures.

Figure 3 below illustrates the procedure for mechanical automated disaggregation using a commercially-available disaggregation machine.


Figure 3 (from Reference 2). Automated method. (A) Slicing tissue into smaller fragments, (B) Placing tissue into commercially available disaggregation device, (C) Device in disaggregation machine, and (D) Cell suspension aspirated into tuberculin syringe.

Study results are mixed on the utility and yield of these types of automated mechanical disaggregation techniques. Novotny et al.⁷ used this technique on non-aspirable bone marrow samples and found at least equivalent yield compared to the vortex technique described by Vos et al.⁶ However, other studies have demonstrated the limitations of the automatic mechanical disaggregation method (see histograms below in Figure 4) and obtained higher percentages of diagnostic flow samples utilizing manual techniques, particularly in large cell lymphomas.^{2, 6}


Figure 4 (from Reference 2). Comparison of flow histograms obtained using manual versus automatic tissue disaggregation.

Table 1. Summary Table of Advantages and Disadvantages.

Technique	Advantages	Disadvantages
Needle, Mesh, Scraping, Mincing, Tissue press	1) Gentler tissue processing 2) Higher cell yield 3) Decreased cell debris	1) Increased tissue handling 2) Limited by necrosis and fibrosis
Vortex	1) Allows histologic and flow cytometric correlation on same tissue sample 2) Well-suited for extremely small samples 3) Minimal tissue handling 4) For use with tissue and non-aspirable bone marrow biopsies	1) Limited by necrosis and fibrosis
Medimachine	1) Commercially available 2) Standardized 3) Operator-independent 4) Closed disposable system	1) High shear forces 2) Particularly susceptible to loss of large cell lymphomas

VIABILITY

Assessment of cellular viability is essential for accurate flow cytometric analysis and may be a good way to demonstrate the effectiveness of your preferred tissue disaggregation technique. Non-viable cells can demonstrate non-specific antibody binding of fluorescent antibodies and/or display unusual autofluorescence that can interfere with identification of abnormal populations. In addition, in quantitative assays such as lymphocyte enumeration (T-cell subsets) or CD34+ stem cell counts, the presence of non-viable cells can cause inaccurate results. Thus, the purpose of viability assessment is to ensure sample integrity, allow for precise measurement of antigenically-intact (viable) cells, and quantitation and/or exclusion of non-viable cells. In addition, CAP accreditation requires “a policy for determining when the percentage of viable cells in each test specimen should be measured” (refer to the CAP Flow Cytometry Checklist; FLO.30610).⁸⁻⁹

Cellular viability is often also closely related to specimen type. For example, fresh peripheral blood samples normally demonstrate excellent viability, which decreases over time. However, decreased cell survival is common in tissue samples and increases in high-grade neoplasms with high cellular turnover or necrosis or those exposed to tissue processing or prolonged storage.

There are several methods for assessment of cell viability. The two most commonly used fluorescent viability dyes are propidium iodide (PI) and 7-aminoactinomycin D (7-AAD)¹⁰. Both of these dyes intercalate with (attach or insert into) DNA base pairs meaning that cellular staining (positivity) is indicative of a compromised cellular membrane. Both dyes are excited by a blue

laser (488 nm) with emission maxima of 617 nm (PI) and 647 nm (7-AAD).⁸⁻¹⁰ Gating on viable events is performed by drawing a gate around the PI or 7-AAD negative events on viability dye versus side scatter dot plots. A rough estimate of viable cells can also be approximated by excluding low forward scatter (FSC) and high side scatter (SSC) events (see Figure 5 below).¹⁰


Fig. 1: Viability gating with 7-AAD


Fig. 2: Non-viable cell exclusion by FSC/SSC

Figure 5 (from reference 9). Viability gating

SUMMARY

Tissue disaggregation remains one of the most important and problematic steps in solid tumor analysis by flow cytometric immunophenotyping. The method used for tissue disaggregation can have a definite impact on cell loss, viability, non-specific antibody binding, and antigenicity. There are several options for performing tissue disaggregation available as discussed above. The goal here is not to present a single optimal approach, but to provide information regarding the advantages and disadvantages of each technique.

References:

1. Bertram, HC, Check IJ, & Milano, MA. Immunophenotyping large B-cell lymphomas. Flow cytometric pitfalls and pathologic correlation. *Am J Clin Pathol* 2001; 116:191-203.
2. Vallangeon, BD, Tyer, C, & Lagoo, A. Improved Detection of Diffuse Large B-Cell Lymphoma by Flow Cytometric Immunophenotyping – Effect of Tissue Disaggregation Method. *Cytometry B Clin Cytom*, 2015 Sept 9; 90(5):455-61.
3. ICCS Education committee/Q&A by G Deeb and A Illingworth "What are the optimal methods for processing tissues for flow cytometry?"
4. Cornacchiari, A, Grigolato, PG, Facchetti, F, Morassi, ML, Cadei, M, Alpi F, Battocchio, S, & Chirico, E. Usefulness of the scraping method for DNA flow cytometry in breast tumors. *Cytometry*. 1995 Mar 1; 19(3): 263-6.
5. Frank X. Torres, Patricia G. Mackowiak, Ron D. Brown, Michael D. Linden, Richard J. Zarbo, Comparison of Two Methods of Mechanical Disaggregation of Scirrhous Breast Adenocarcinomas for DNA Flow Cytometric Analysis of Whole Cells, *American Journal of Clinical Pathology*, Volume 103, Issue 1, 1 January 1995, Pages 8–13.

6. Vos, J, Simurdak, J, Davis, B, Myers, J & Brissette, M. Vortex disaggregation for flow cytometry allows direct histologic correlation: A novel approach for small biopsies and inaspirable bone marrows. *Cytometry B Clin Cytom*, 2003, 52B:20-31.
7. Novotny JR, Schmücker U, Staats B, Dührsen U. Failed or inadequate bone marrow aspiration: a fast, simple and cost-effective method to produce a cell suspension from a core biopsy specimen. *Clin Lab Haematol*. 2005 Feb;27(1):33-40.
8. ICCS Ask an Expert by Sarah Sivers & George Deeb. "How do I establish and validate a panel using a viability dye such as 7-AAD? How do I perform compensation for this panel?"
9. ICCS Ask an Expert by Adam Seegmiller "How and in what situations should viability of specimens for flow cytometry be assessed?"
10. Schmid I, Krall WJ, Uittenbogaart CH, et al. (1992) Dead cell discrimination with 7-amino-actinomycin D in combination with dual color immunofluorescence in single laser flow cytometry. *Cytometry*. 13:204-208

For any questions on this module or any other suggestions, please email info@cytometry.org

Reviewed and approved by: George Deeb, MD and Jean Oak, MD

The documents posted on ICCS website may contain product or vendor names which are provided for platform specific guidance. Any reference within the ICCS Quality and Standards modules to any vendor, product or educational material by trade name, trademark or manufacturer does not constitute or imply the endorsement or recommendation by ICCS.