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Quantitative Analysis of Phospholipids by HPLC with a Light Scattering Evaporating Detector –

Summary

A method is described for separation and quantitative analysis of various phospholipids and sphingolipids by high performance liquid chromatography on a silica gel column with a mobile phase consisting of chloroform, methanol, and ammonium hydroxide. Gradient elution is necessary.

Key Words:

High performance liquid chromatography
Phospholipids
Light scattering evaporative detector

2.1 Equipment

The chromatographic system comprises the following components:

- 2 pumps LC6A (Shimadzu)
- 1 gradient controller SCL6A (Shimadzu)
- 1 650 μ l mixing chamber for pumps LC6A (Shimadzu)
- 1 autosampler SEDEX 100 (SEDERE, 8 Rue Eugène Henaff, 94403 VITRY-SUR-SEINE Cedex B.P.52, France)
- 1 light-scattering evaporative detector (SEDERE)
- 1 "Autofix II" system fitted with a 125 mm "Lichrocart" cartridge filled with "Lichrospher" Si 60 (5 microns), fitted with a guard column with the same material (Merck)


The chromatographic system is equipped with an on-line filter (Waters) without dead volume between the injector and the column. The system allows successive injection of many samples according to a perfectly standardized cycle.

2.2 Reagents and Standards

- chloroform p.a. (Carlo Erba)
- methanol p.a. (Carlo Erba)
- ammonium hydroxide at 30 % (Carlo Erba)

Phospholipids and cerebroside standards were purchased from Sigma.

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

2.3 Elution Gradient

000121

Separation of a standard mixture of nine (phospho) lipids. 1, CR; 2, CROH; 3, PE; 4, PI; 5, PS; 6, PC; 7, PA; 8, Sph double peak; 9, LPC. Conditions in text.


Figure 3. Calibration curves for lipid standards

	AIRES							
	injection 1	injection 2	injection 3	injection 4	injection 5			
PIC A	$2.09 \cdot 10^9$	$2.12 \cdot 10^9$	$2.00 \cdot 10^9$	$1.88 \cdot 10^9$	$1.91 \cdot 10^9$			
PIC B	$0.47 \cdot 10^9$	$0.50 \cdot 10^9$	$0.53 \cdot 10^9$	$0.49 \cdot 10^9$	$0.51 \cdot 10^9$			
PIC C	$1.94 \cdot 10^9$	$1.87 \cdot 10^9$	$2.00 \cdot 10^9$	$2.13 \cdot 10^9$	$2.06 \cdot 10^9$			

000123